BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 0719004306	Date: 04DEC24	Proofread No.: 05
Designer: Z.B.	Page: 1 of 2	
Colour Reference: BLACK		

Baxter

REF RMC4916

Irrigation jet (en)

(1) Protective Cap (2) Spike

Approx. Length: 19 cm

Intended Purpose: For the delivery of irrigation solutions from the fluid container to the irrigation site.

Intended Users: Healthcare Professionals

Sterile and Non-pyrogenic

Directions for Use: Use Aseptic Technique

1. Remove Protective Cap (1) from Spike (2).

2. Insert spike (2) into irrigation solution bag.


3. Fully prime the set to remove all air bubbles from set prior to use.

Cautions:

- Do not use if package has been opened or damaged or if tip protectors are loose or missing.
- The set can be used under pressure to assist flow (max. 300mmHg).
- Do not allow air to be trapped in set.
- Do not remove from pouch until ready to use.
- Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.
- Single use only. Do not resterilize.
- Not I.V. compatible.

Notes:

- All incidents should be reported to the manufacturer as identified on the label. In case of incidents involving death or serious injury report them to the competent authority as well.
- Store at room temperature.
- Replace set every 72 hours or as per institutional protocol, whichever comes first.
- Dispose as per healthcare provider's policy.
- The fluid path of this product contains ABS and PVC. Ensure that the drug is compatible with these materials.
- This product is not manufactured with natural rubber latex.
- Contains less than 0.1% w/w DEHP
- See the glossary for all symbols and definitions.

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Baxter is a registered trademark of Baxter International Inc.
Baxter è un marchio registrato di Baxter International Inc.